
















Available Safe Patient Handling (SPH) Equipment


How to Access: To provide patient comfort and prevent employee injuries, SBM has available the following safe patient handling equipment. Equipment along with operating instructions is obtained through Capacity Management (and delivered by Distribution Services Transport staff) or from Lawson, depending on the type of equipment. Contact your Clinical Educator for training. Refer to the manufacturer's website for additional information.




Equipment/Quantity	Manufacturer/ Model #	Capacity (lbs.)	How to Request/Order	Photograph	Manufacturer's Website (right click on website link; select "Open Hyperlink")
General SPH/Bariatric Equipment					
Compella Bariatric Bed (6)	Hill-Rom P7800	1,000	Capacity Management		http://www.hillrom.com/usa/Products/Category/Hospital-Beds/compella-bariatric-bed/
TotalCare Bariatric SpO2RT Bed (4)	Hill-Rom P1840DCA0000 10	500	Capacity Management		https://www.hillrom.com/en-us/products/smart-beds-and-surfaces/totalcare-spo2rt-2-icu-bed/
Progressa Bed System (7)	Hill-Rom P7500A000319	500	Capacity Management		https://direct.hillrom.com/hillromUS/en/Products/Beds/PROGRESSA-BED/p/388B186E-164E-4384-8E1A-B9BFC3AC6B3
Stretcher Chair (Cardiac Chair) (10)	Stryker Model #550	400	Capacity Management		https://www.stryker.com/us/en/acute-care/products/stretcher-chairs.html






Equipment/Quantity	Manufacturer/ Model #	Capacity (lbs.)	How to Request/Order	Photograph	Manufacturer's Website (right click on website link; select "Open Hyperlink")
General SPH/Bariatric Equipment (continued)					
Bariatric Wheelchair (3)	Hill-Rom P654A SD100370	800	Capacity Management		https://hillrom.com/en-us/products/
Standard and Bariatric Wheelchair (~100 standard and 6 bariatric)	MDS806950 MDS809950 Medline	500 850	Capacity Management		https://www.medline.com/product/Extra-Wide-Wheelchairs/Standard-Wheelchairs/Z05-PF05040?question=MDS806950&index=P1&indexCount=1 https://www.medline.com/product/Shuttle-Extra-Wide-Wheelchairs/Standard-Wheelchairs/Z05-PF05041?question=MDS809950&index=P1&indexCount=1#
Bariatric Shower Chair with Backrest (6)	Medline G2-102BX1	550	Capacity Management		https://www.medline.com/product/Shower-Chair-with-Back/Benches/Z05-PF04624?question=&index=P9&indexCount=9
Bariatric Commode Chair (3)	Drive 11117N-2	650	Capacity Management		http://www.drivemedical.com/index.php/bariatric-folding-commode-18.html
Bariatric Commode Chair (4)	Medline Guardian MDS89668XW	850	Capacity Management		http://www.medline.com/product/Bariatric-Drop-Arm-Commode/Commodes/Z05-PF04649?question=MDS89668XW&index=P1&indexCount=1
Bariatric Commode Chair (2)	Hill-Rom SC100398	1000	Capacity Management		https://hillrom.com/en-us/products/

Equipment/Quantity	Manufacturer/ Model #	Capacity (lbs.)	How to Request/Order	Photograph	Manufacturer's Website (right click on website link; select "Open Hyperlink")
Folding Walker with wheels (~60)	MDS86410W54B MDS86410XWWB Medline	400 (~30) 600 (~30)	Capacity Management		https://www.medline.com/product/Two-Button-Folding-Walkers-with-5-Wheels/Walkers/Z05-PF05011?question=MDS86410W54B&index=P1&indexCount=1 https://www.medline.com/product/Bariatric-Folding-Walkers/Walkers/Z05-PF05016?question=MDS86410XWWB&index=P1&indexCount=1
Bariatric Scale (1)	Scale-Tronix (Welch Allyn) 6006	600	Capacity Management		https://www.welchallyn.com/en/products/categories/physical-exam/scales/portable/scale-tronix-portable-scales.html
Bariatric Wheelchair Scale (1)	Detecto 6550	1000	Capacity Management		https://www.detecto.com/product/product-overview/WheelChair-Scales/6550-Portable
Bariatric Gurney/Autopsy (1)	Ferno 24 MAXX 236- 1417-OOC	1000	Capacity Management		http://www.fernomortuary.com/en/Products/First%20Call/cots/24-maxx-mortuary-cot.aspx?ec_trk=followlist&ec_trk_data=cots
Bariatric Cadaver Carrier (2)	Mortech (Medline) MDR6000039H	600	Capacity Management		http://www.medline.com/sku/item/MDPMDR600039H;ecomsessionid=xTAAMGhn5bsEiBc-IVY3+Q?skuIndex=S4&question=&flowType=&indexCount=
Sit to stand					
Sara Steady (25)	ArjoHuntleigh NTB2000	400 Client Height 4'11' – 6'6"	Capacity Management		https://www.arjo.com/int/products/patient-handling/standing-and-raising-aid/sara-steady/






Equipment/Quantity	Manufacturer/ Model #	Capacity (lbs.)	How to Request/Order	Photograph	Manufacturer's Website (right click on website link; select "Open Hyperlink")
1- Lotus XL Active stand – up aid	Lopital SKU 59008270	500 Client Height 5'7" – 6'5"	Distribution		https://www.lopital.com/products/manual-transfer-aids/lotus-xl-active-stand-up-aid/ 
TruRize Powered Clinical Chairs (6) 1– Pavilion ICUs 1– 9 South 1– 17 South 1– 18 South 1-13N 1-19N and 19S Share 1-9W	Stryker 3333-000-000	350	Capacity Management		https://patientcare.stryker.com/en/products/clinical-seating/clinical-seating/trurize
Electric Stand Assist Lift (1)	Medline MDS400SA	400	Capacity Management		http://www.medline.com/product/Electric-Stand-Assist-Lift/Lifts/Z05-PF04443?question=mds400sa&index=P1&indexCount=1


Equipment/Quantity	Manufacturer/ Model #	Capacity (lbs.)	How to Request/Order	Photograph	Manufacturer's Website (right click on website link; select "Open Hyperlink")
Electric Stand Assist Lift (1)	Medline MDS600SA	600	Capacity Management		http://www.medline.com/product/Electric-Stand-Assist-Lift/Lifts/Z05-PF04443?question=mds600sa&index=P1&indexCount=1

Lateral Transfer and/or Repositioning or Lift from Floor					
Equipment/Quantity	Manufacturer/ Model #	Capacity (lbs.)	How to Request/Order	Photograph	Manufacturer's Website (right click on website link; select "Open Hyperlink")
Hovermatt Air Transfer System (1)	Hovertech International HM34HS	1,200	Capacity Management		http://www.hovermatt.com/hovermatt.html
E-Z Mattress (2)	E-Z Way, Inc. Regular	1,200	Capacity Management		https://www.ezlifts.com/lateral-transfer-solutions
Antistat Patient Shifter slide board (~20) Note: for lateral transfer Only	Alimed Inc. NA	250	Typically stored in clean utility rooms		http://www.alimed.com/patient-shifters-antistat-and-conductive-12239.html

HoverJack Air Patient Lift (1) Note: for lifting patient from floor	Hovertech International HJ3902T	1,200	Call Distribution Services at 4- 2980		http://www.hovermatt.com/hoverjack.html
ComfortGlide TRP Air Assisted Repositioning Sheets	Medline Comfort Glide Air MSC60011	1000lbs	Storeroom Lawson #194288 Deposit in marked bag soiled utility room.		https://www.medline.com/skin-health/comfortglide/#single/0 
Medline Wide Transfer Belt with handles	Medline MDT821203BL CH	NA	Capacity Management		https://www.medline.com/product/Medline-Wide-Transfer-Belt-with-Handles/Belts/Z05-PF80076
Equipment/Quantity	Manufacturer/ Model #	Capacity (lbs)	How to Request/Order	Photograph	Manufacturer's Website (right click on website link; select "Open Hyperlink")
Repositioning					
Medline ComfortGlide TRP Air Assisted Repositioning Sheets	MSC60011	1000LB	Capacity Management		<p>SHEET SYMBOL IDENTIFICATION</p>  <p>1. Flap over compartments 2. Lamin bed handles 3. Use control rod and head handles to center patient.</p> <p>4. Roll and secure on floor using roller sheet.</p> <p>5. Reposition patient with padding down.</p> <p>Air Blower Nozzle Indicator Shows direction to point blower over ComfortGlide® Air Sheet.</p> <p>Head Placement with ComfortGlide Indicates where the head is placed at the top of the Air Sheet.</p> <p>WARNING The Medline ComfortGlide® Air Sheet has been determined through scientific studies to be safe. Safe use requires this product and accompanying ComfortGlide® Drypad be constructed of the correct size and used for the intended purpose of lifting, moving, and repositioning. An abnormal sheet product may lead to the accident or injury during the lifting. The ComfortGlide® Air Sheet is not to be used in MRI safe zone in available open magnet. Please follow your facility protocol for using the Blower in an MRI environment. For more details, review related ComfortGlide® Drypad/Sheet imaging.</p> <p>CLEANING PROCEDURE The ComfortGlide® Air Sheet is single patient use only. Do not handle the ComfortGlide® Air Sheet or ComfortGlide® Drypad. The ComfortGlide® Drypad should be replaced and discarded when called to replace the ComfortGlide® Air Sheet. Only approved personnel should use the sheet. If the ComfortGlide® Air Sheet becomes discolored or heavily soiled, discard sheet. Follow all infection control guidelines.</p> <p>RETURNS For returns and concerns, please contact your Medline representative.</p>  <p><small>© 2019 Medline Medline is a registered trademark of Medline Industries, Inc. All other trademarks are the property of their respective owners.</small></p>
Medline Comfort Glide Repositioning Sheet with	Medline CGLIDESHEET CGSHEETXLD	575 (std) 1,000 (XL)	Lawson in- house 138311 (std)		http://www.medline.com/product/Comfort-Glide-Repositioning-Sheet/Drawsheets-and-Positioning-Sheets/Z05-

handles and drypad, no wedges (single patient use)			138309 (XL) 138308 (std pad) 138313 (XL pad)		PF80216?question=cglidesheet&index=P1&indexCount=1
Single Patient SliderZ perforated transfer sheet (single use)	Medline MDS30180	400	Lawson in-house 109988		https://www.medline.com/product/Disposable-Transfer-Sheet/Drawsheets-and-Positioning-Sheets/Z05-PF113929;ecomsessionid=SC7BemulCeCUEyHPms1pVw?question=transfer+sheet&index=P1&indexCount=1
Ceiling Mounted Lifts					
Bariatric & Metabolic Weight Loss Center, 9N-010 (1)	ArjoHuntleigh Maxi Sky 600	600	NA		https://www.arjo.com/en-us/products/patient-handling/ceiling-lift/
Hospital Pavilion 2 lifts per in-patient floor, 6W-9W Rooms 0802 and 0201 (8)	Handicare (Prism Medical) C625	625	NA		https://www.handicareusa.com/product/c-625-fixed-ceiling-lift/
Autopsy cadaver lift 2-736 (1)	Guldmann GH3 350 120 1110	550	NA		https://www.guldmann.com/us/products/ceiling-lift-systems/lifting-motors/gh3

Equipment/Quantity	Manufacturer/ Model #	Capacity (lbs)	How to Request/Order	Photograph	Manufacturer's Website (right click on website link; select "Open Hyperlink")
Mobile Lifts & Slings					
Golvo Patient Lift (16)	Hill-Rom 8008	440	Capacity Management		https://www.hill-rom.com/globalassets/website-documentation/english-websites-us-int/patient-handling-us-int/mobile-lifts/golvo-80008008-us/golvo_instruction_guide_7en140107-03.pdf
Liko Viking XL Patient Lift (6)	Hill Rom	660	Capacity Management		https://www.hill-rom.com/international/Products/Products-by-Category/patient-handling/mobile-patient-lifts/viking-xl-patient-lift/
Medline Electric Bariatric Lift (6)	Medline MDS600EL	600	Capacity Management		http://medline.com/irj/catalog/details/MD_S600EL/
E-Z Way Smart Lift (2)	E-Z Way, Inc. 598	600	Capacity Management		https://www.ezlifts.com/mobile-lift-solutions
E-Z Way Smart Lift (2)	E-Z Way, Inc. 498	1,000	Capacity Management		https://www.ezlifts.com/mobile-lift-solutions

Equipment/Quantity	Manufacturer/ Model #	Capacity (lbs)	How to Request/Order	Photograph	Manufacturer's Website (right click on website link; select "Open Hyperlink")
Mobile Lifts & Slings (continued)					
Various Slings (e.g., 4 points sling)	Medline, EZ Lift, Golvo	Refer to specific lift	Capacity Management Return in bag to Linen for laundrying	Example of a sling: 	Hill-Rom information: http://www.hill-rom.com/globalassets/website-documentation/english-websites-us-int/patient-handling-us-int/slings-and-lift-sheets/soft-originalsling-us/basic_models_quick_reference_guide_7_en161501-02.pdf

Revision 4/2026